



	FAB SWAB	Intergriswab	Cap-Shure
<p>Protects With Hard Shell i.e. Has a hard shell container that closes tightly and protects both the entire shaft and the cotton tip of the swab</p>	Yes	No	No
<p>Air Dries Through Breathable Filter i.e. air dries through a Tyvek breathable filter that covers 8 large integrated air holes</p>	Yes	No	No
<p>Avoids contamination with the hands of examiner i.e. has a specially designed handle so that the swab can be handled without touching the shaft</p>	Yes	No	No
<p>Avoids contamination after sampling i.e. no unprotected air holes to contaminate or cross contaminate of the scene, other evidence, equipment or the laboratory.</p>	Yes	No	No

Handle your most valuable evidence and sensitive traces with FAB-SWAB™