



# FAB SWAB™

Filtered Air Breathable Swab

**A swab that protects, air dries, avoids contamination, does not require drying equipment**

## Protects

Has a hard shell container that closes tightly, protects entire shaft and the cotton tip of the swab.

## Air Dries

Air dries it with built in Tyvek® breathable filter covering 8 large integrated air holes.

## Optimizes Drying

The design, size and the placement of the filter protected air holes optimizes drying.

## Avoids Contamination

No unprotected air holes such as those found in Integriswab, Cap-Shure swab.

**Visit [fabswab.com](http://fabswab.com) for more details**



Abacus Diagnostics • 6520 Platt Ave #220 • West Hills, CA 91307 • USA  
 Phone (818) 716-4735 • Toll Free (877) 225-9900 • Fax (818) 716-9471  
 Email: [CustomerService@abacusdiagnostics.com](mailto:CustomerService@abacusdiagnostics.com) • [www.abacusdiagnostics.com](http://www.abacusdiagnostics.com)